


Automation - Integration

Last Modified on 2021-05-20 18:27

To automatically send examinations to an EMR, select an option under **Integrations** has the Automation Type when [creating your automation rule](#):



Select Automation Type

Athena Health ▼

Integrations

- Allscripts Professional EHR Integration
- Allscripts Touchworks EHR Integration
- Athena Health
- REDCap


After the process has been created, the configuration screen displays.

First, select **Add New Condition** to dictate what triggers the automated process. It is important to note that a condition is not required for creating an automatic process. **If conditions are not added, the automatic process is always applied.**

Select the red X to remove a condition.

If adding a condition, select an option from the **When** drop-down menu to specify when the automation is applied:

When should this automation occur?



Add new condition

When Status ▼ Equals [] X

Study

- Status

Instance

- Destination AE Title
- Destination AE Title Extension
- Is or contains a PDF

DICOM Fields

- Choosing an option under **Instance** runs the process for specific instances
 - **Destination AE Title** or **Extension** will run the automated process when a specific imaging system is used
 - Type the name of the imaging system into the Equals text box
 - For example, if a specific machine is used for ultrasounds, type the name of the machine into the text box. Each time the machine is used, the examination will be sent to the patient's chart within your EMR
 - **Is or contains a PDF** will run the automated process when the study contains a PDF document
 - Leave the Equals text box empty
- Choosing an option under **DICOM Fields** runs the process when specific information is entered into that DICOM field
 - Type the information that needs to be entered into the Equals text box
 - For example, if Dr. Julie does all patient exams, choose Operators Name as



- For example, if Dr. Julie does all patient exams, choose Operator's Name as the DICOM field and type Dr. Julie into the text box. Each time Dr. Julie is entered into the Operators Name field of the imaging machine, the examination will be sent to the patient's chart within your EMR

Currently, Tricefy works with Allscripts, athenahealth, and Redcap.

- Select **Automatically Generate Outbound HL7 Message** to integrate with any EMR using HL7.

Select the appropriate article below to finish setting up your automation rule:

 [athenahealth users](#)

 [HL7 message configuration](#)
