




Tricefy Uplink Requirements

Last Modified on 2022-04-07 10:04

Tricefy Uplink has minimal requirements:

1. An account with [Tricefy](#) (user name and password)
2. DICOM-capable and DICOM-enabled Ultrasound machine(s)
3. Internet connection with “broadband speed” or better.
 - **Tip:** Faster internet connections will result in faster uploads
4. DICOM -licensed ultrasound system(s) connected to an active network
 - The ultrasound system and Uplink software must be connected to the same network
5. A computer, server, or virtual machine for running the software. Tricefy Uplink can be run in conjunction with other applications and does not need a dedicated computer. A static/fixed IP address is strongly recommended.
 - a. Window Requirements:
 - i. Windows 10 (previous versions are not HIPAA compliant) with a 2012+ server **
 - ii. 40MB of disk space necessary for processes and data files
 - iii. Computer with two or more CPU cores (dual core or better)
 - iv. [Static IP Address](#)
 - b. Mac Requirements:
 - i. OS X 10.6
 - ii. 40 MB of disk space necessary for processes and data files
 - iii. Computer with two or more CPU cores (dual core or better)
 - iv. [Static IP Address](#)
5. Up-to-date modern browsers for viewing Tricefy on a computer or mobile device
 - a. Supported browsers
 - Chrome
 - Firefox
 - Edge
 - b. Using an out-of-date browser will redirect to the mobile version of Tricefy for security purposes

Note: Internet Explorer will no longer be supported after April 2022

 Learn more: tricefy.help/help/ie

** If you are unsure whether your Windows Server is 2008 or 2012, refer to your password. 2012 servers require a 12-character minimum password with complex password rules.
