

Last Modified on 2024-09-24 15:55

To create a study with DICOM, upload the examination by selecting Upload DICOM from + menu at the top of the Study List or Consults screen:



If you are creating a study without DICOM, such as a study that only contains PDF documents, refer to the link below:

Creating a Study without DICOM

When uploading DICOM, a pop-up window displays so you can select files from your computer/network:

Back Upload DICOM Files									
Account BN T	inprogress v Filename	Туре	Size	Status	DICOM	Progress			
Files Choose Files No file chosen									
Begin Upload									

• If you are a member of multiple clinics, confirm the clinic shown in the Account field is correct. To change the clinic, choose another option from the drop-down menu.

Next, select Choose Files to browse to a location on your computer/network and select your files. Selected files will display under File Name.

inprogress •					
Filename	Туре	Size	Status	DICOM	Progress
20160825-KRISTI.dcm		290760	pending	Yes	

- The Status column indicates if the file is Uploaded, Pending upload (either because the Begin Upload button was not selected or because the file is in queue) or Failed to upload.
- The Progress column displays a progress bar that turns blue as the file is uploaded
- Note: single file size limit is 300 MB

Tricefy will create the study from the file(s) and use the date, patient name, and other information that was entered into the imaging system when performing the exam.

Verification

To make sure the study was created, check the My Activity tab of the side information panel on your study list. It should say a study was created and provide a link.



