


Step 1: Activate your Voluson


Last Modified on 05/04/2017 5:14 pm EDT

The following steps will have your Tricefy activated in just seconds!

- 1) Press the **Utilities** button on the front panel (this is a physical button), and select **Setup** on the touch screen:




- 2) Select **Connectivity**:



- 3) Select the **Tricefy** tab and activate your account.

Note: The Tricefy tab is only available on new machines that have Tricefy Inside.





A - Enter your email address to authorize the account. This should be an IT or clinic administrator.

B - Select **Activate Account**. This process may take up to 20 seconds.

- A **green** message means the activation was successful
- An **orange** message the activation failed
 - If activation fails, confirm there are not any firewalls or routers blocking Port 443
 - If activation continues to fail, select the **Contact Tricefy Help Center** button at the bottom of this screen and we will contact you shortly to assist in getting your account activated

C - Select **Test Connection** to confirm your ultrasound system can connect to Tricefy

- A **green** message means the test was successful and you can continue with the steps below
- An **orange** message the ultrasound system failed to connect to Tricefy
 - If the test fails, confirm there are not any firewalls or routers blocking Port 443
 - If the test continues to fail, select the **Contact Tricefy Help Center** button at the bottom of this screen and we will contact you shortly to troubleshoot the issue and get your system connected.

If you have multiple ultrasound systems, each system requires activation (and all will link to your Clinic's Tricefy account).

 [Move to Step 2: Configuration](#)
