Tracefy Correcting a Patient ID

Last Modified on 2025-03-04 12:19

Note: If the patient is not already located in Tricefy (this is their first exam), you can fix the patient ID by editing the patient information instead.

Instructions for editing patient information (adding a unique ID for a new patient)

If exams have been uploaded for the patient in the past, continue with the instructions below for correcting their patient ID.

Studies uploaded to Tricefy using the wrong patient ID can be transferred to the correct patient ID.

1) Select the : (action menu) at the top of the screen and choose Edit Study:

Studies Open In Viewer									
		Full Name	Study Date ~	Status	Study Details	Interface Status		Anonymo	Attach File(s) Edit Study
Ŧ	•	Dean Gerald	2018-02-21 5:21 PM MST	None		No status	- P	IFO	↔ Move
T	•	Anonymous	2018-02-21 4:14 PM MST	None	I 1	Fail		Study Da 2018-02-2	Download Studies Delete
T	•	Anonymous	2018-02-21 1:05 PM MST	Ongoing	⊞ 2 9 1	Success	Í	Study Ex 2043-02-2	Clear Selection

2) Search for the correct patient (using their name or correct patient ID) and select the appropriate result:

I	Edit Study	×					
Patient							
	Thurman Gerald Dean (fd0d0a29-5f2a-4c43-aed1-061ce4a7e36a)						
	0045						
	BIELENBERG DAVID (004540274)						

Note: If the patient is not already located in Tricefy (this is their first exam), fix the patient ID byditing the patient information

After selecting the correct patient, you will notice the study is now assigned to the right patient ID.

The information that was uploaded to Tricefy (the data in your ultrasound system) is not changed. This data can only be corrected by editing the information using your ultrasound machine.

However, the corrected changes are reflected when downloading the study from Tricefy.

Instructions for downloading content from Tricefy