



Additional Uplink Documentation (technical)

Last Modified on 2023-08-22 12:17

You can change the default behavior of Uplink by entering an expert code during installation:

| Expert Setting | Behavior |
|----------------|---|
| CP | Cloud processing |
| Cached | No local processing |
| LPU-O | Local processing, optimized pixels - Default for trice4 * |
| LPU-D | Local processing, deflated pixels |
| LPU-R | Local processing, retain original pixel data - (former triceU) ** |
| Anon-LPU-O | Anonymized strict, optimized pixels - (Former triceN)*** |
| Anon-LPU-D | Anonymized strict, deflated pixels |
| Anon-LPU-R | Anonymized strict, retain retain original pixel data |
| AnonNS-LPU-O | Anonymized non-strict, optimized pixels |
| AnonNS-LPU-D | Anonymized non-strict, deflated pixels |
| AnonNS-LPU-R | Anonymized non-strict, retain retain original pixel data |
| Pseud-LPU-O | Pseudonymization strict, optimized pixels |
| Pseud-LPU-D | Pseudonymization strict, deflated pixels |
| Pseud-LPU-R | Pseudonymization strict, retain original pixel data |
| 4DVIEW | Setting needed when using 4DVIEW |

Anonymized data cannot be re-associated with the original PHI, whereas pseudonymization can be tracked back. This application is most commonly used in clinical studies.

 [Learn about strict vs non-strict anonymization](#)

To apply an expert setting, you must install the [Advanced Uplink](#).

File Delay

Enabling file delay will delay sending any content to Tricefy until the examination is complete. To activate this feature, install the [Advanced Uplink](#).

Select File Export & DICOM on (third screen of the installation wizard) and enter one of the following two options into the File Extension text box:

.* Delay : Entering this option will delay sending any file type to Tricefy

.PDF Delay : Entering this option will delay sending only PDFs to Tricefy (exam images will be sent as normal)

Note: Make sure there is a space between the suffix (*) and "Delay" or "PDF Delay"
