



Templates - Adding and removing Headers

Last Modified on 08/15/2019 4:18 pm EDT

To get started customizing your template, it is recommended to begin with the Tricefy global template that best fits your needs. Changes to your template are automatically saved. The template editor should not be used with tablets and mobile devices.

Titles and headers are used to organize the templates. Titles typically go at the top of the page and describe the report; headers are typically used at the start of a new section to describe that section.

Adding a Title

To add a title:

- 1) Hover your mouse where you want the title (typically at the top of the page)
- 2) Select the + button to add a new item
- 3) Choose **Text** from the menu, followed by **Title**

Adding a New Header

- 1) Hover your mouse where you want the header
- 2) Select the + button to add a new item
- 3) Choose **Text** from the menu, followed either by **Header** Or **Header with Underline**

Once you've added a new title or header, you can adjust the color, size, and alignment using the Text Settings in the right sidebar:

The screenshot displays the Tricefy template editor interface. On the left, a 'Biophysical Profile' section is visible, containing a table with fetal assessment data and a text box for 'Fetal Anatomy and Radiometry'. The right sidebar, titled 'TEXT', provides settings for the selected text element, including font size (18), bold (B), italic (I), and background color. Below the 'TEXT' section, the 'ELEMENT' section shows a 'Background' setting with a color picker.

Biophysical Profile				
Qualitative AFV	2	Growth Body Movement	2	Total 8
Fetal Tone	2	Fetal Breathing Movement	2	Empty

Fetal Anatomy and Radiometry

Fetal anatomic and biometric evaluations were not performed on today's exam.

Text

TEXT

18 B I

ELEMENT

Background

Removing a Header

To remove a header,

- 1) Click on the Header to help locate it on the left sidebar
- 2) Hover your mouse over the header on the left sidebar
- 3) Click the trashcan icon that appears when your mouse is hovered:

