



## Exporting all data out of Tricefy

Last Modified on 2026-06-17 19:14

You can export all data for storage to another system.

- 1) Contact your sales representative or our support center ( [support@triceimaging.com](mailto:support@triceimaging.com)) to start the process
  - 2) The Trice team will assist installing the [Tricefy Downlink](#) software on a computer in the customer network and help configure the directory to write the data into
  - 3) The Trice team will then queue up all study data in the Tricefy cloud account(s) for transfer through to the local downlink
  - 4) Once all export jobs have completed, the customer and the Trice team will verify all data has been transferred (usually by confirming study counts between systems)
-