


## Tricefy Uplink (Installation and more)

Last Modified on 2021-03-18 18:31



If you are using a new GE Voluson, you may have [Tricefy Inside](#). This feature is a pre-installed Uplink - all you have to do is activate your Tricefy account.

Look for the [Tricefy Inside](#) sticker on your machine.

## Downloading & Installing Tricefy Uplink

The Uplink most suitable for your clinic was determined when you signed up for your account:

### Typical Uplink

Suitable for most small and medium sized clinics - DICOM is anonymized before it's sent to the patient. It is simple and user friendly; no IT resources are needed

[👉 Typical Uplink Installation Instructions](#)

### Advanced Uplink

Typical Uplink with additional features for integrating with a PACS or EMR:






- Sending studies to Tricefy using ViewPoint
- Routing studies from Tricefy to an external server
- Query and retrieve using a local server

[👉 Advanced Uplink Instructions](#)

If you are performing clinical studies or have a need for additional anonymization, contact our Customer Success Team (by selecting the [Contact](#) button below) to inquire about our Anonymized Uplink. This Uplink anonymizes data before uploading to Tricefy and discards any non-DICOM content.

## Learn More about Tricefy Uplink:



-  [Tricefy Uplink Requirements](#)
-  [Making sure your IP address is Static](#)
-  [Instructions for a New Computer](#)
-  [Instructions for using Uplink to Download from Tricefy](#)
-  [Learn more about anonymization](#)



Did you ever wonder if you could install Uplink on your Android phone?  
[There is an app for that!](#)

---

Instead of uploading exams to Tricefy, you can also download exams using the Tricefy Downlink

-  [Instructions for installing Downlink](#)