



Tricefy Inside (pre-installed)

Last Modified on 05/04/2017 5:40 pm EDT

New GE Voluson ultrasound systems and new Vscan Extend portable systems have Tricefy pre-installed. To find out if your new system has Tricefy Inside, look for the sticker:



With Tricefy Inside, there is no additional software to download and install - you can begin using your Tricefy trial today!

To get started:

Step 1: Activation

Have your IT department or Network Administrator activate your Tricefy account - activation takes under 30 seconds!

 [See Activation Instructions for the Voluson](#)

 [See Activation Instructions for the Vscan](#)

Step 2: Configuration

Setup your machine to send data to Tricefy either at the end of an exam or with the press of a P button

 [See Configuration Instructions for the Voluson](#)

 [See Configuration Instructions for the Vscan](#)

Step 3: Send your first study

Start a new exam, enter an email address or phone number and share your study instantly!

 [See Instructions for Sending a Study using your Voluson](#)

 [See Instructions for Entering Patient Data in your Vscan](#)
