

Patient App Branding

Last Modified on 02/07/2020 1:05 pm EST

Your patients can view their images with either the [landing page](#) (when using a web browser) or with the Tricefy Mobile app ([available for Android and Apple](#)).

You can customize what your patients see using the Patient App Branding settings ([within your Account Settings](#)):

| Kind | Text or Caption | Target URL | Image | |
|------------|-----------------|-------------------------------|-------|-------------|
| Image | | mailto:mail@yourclinic.com | | Edit Delete |
| Plain Text | Text goes here | | | Edit Delete |
| Button | Press me | https://yourclinicaddress.com | | Edit Delete |

- 1 At any time, select the **Preview** button to see how the mobile app will appear to your patients
- 2 The logo shown on the mobile app is the logo uploaded into your [Image Branding account settings](#)
- 3 Type a heading that will appear at the top of the mobile application
- 4 Select **Add Widget** to add text, images, and buttons:

Add Widget

Widget Type

Plain Text

Button

Image

- Choosing **Plain Text** displays a text box for entering text. Type your desired message, followed by **Save**
- Choosing **Button** displays a window for customizing your button:

Widget Type

Button

Text or Caption

Button/Image Action URL


E.g.: <https://triceimaging.com/> or <mailto:support@triceimaging.com>

Save

- Enter what you'd like your button to say into the **Text or Caption** field, and specify what you'd like your button to do in the **Action URL** field
 - Typing a web address into the **Action URL** field will forward the patient to that website
 - Typing "mailto:" followed by an email address (do not use spaces or quotation marks) will direct the patient to their email service
 - Select **Save** after completing both fields
- Choosing **Image** displays a window for uploading an image - you can choose any image that is saved on your computer or network
 - Similar to the button widget, you can also add an **Action URL** so that image functions as a button
 - Select **Save** after uploading your image and adding an optional action

Reorganizing your App

To change the order of how text, buttons, and images are displayed, select-and-hold a widget while moving it up and down to re-order how they are listed:

| Widgets (drag to reorder) | | | |
|---------------------------|-----------------|---|---|
| Kind | Text or Caption | Target URL | Image |
| Plain Text | Text goes here | | Edit Delete |
| Button | Press me | https://yourclinicaddress.com | Edit Delete |
| Image | | mailto:mail@yourclinic.com |  Edit Delete |

Don't forget, you can add multiples of each kind of widget and select the **Preview** button at any time to see how the mobile app will look.

To remove a widget, select the **Delete** button. To modify a widget, select the **Edit** button.



UPDATE: Apple users require a header, logo, and at least one widget for the patient app landing page to display.
