

Templates - Export Fields

Last Modified on 2020-11-10 17:52

When exporting reports to CSV, specify column names into the **Export Name** field (located at the bottom of the right side panel):

The screenshot shows the 'Export Fields' configuration interface. The central area lists various fields for selection, including Patient Information (Patient ID, Referring Physician, etc.), Fetal Biometry (BPD, HC, AC, FL, etc.), and Amniotic Fluid Index (Q1, Q2, Q3, Q4, AFI). The right-hand panel contains configuration options: DELIMITER (set to '(sentence ending)'), UNIT OF MEASURE (empty), PRECISION (set to 1), DECIMAL SEPARATOR (set to '(dot)'), OBJECT TYPE (set to 'Number'), and EXPORT NAME (a text input field highlighted in red). Below the EXPORT NAME field, there is a note: 'Set this field if you would to include the value when exporting data to integrations.'

If nothing is entered into this field, the column header will either match the mapped field (if possible) or assigned a number (example: "field 10").

Once export fields are entered into a template, they are saved for next time the template is used.

Changing an Export Name

If your export date range includes reports using the template before and after the export name changed, the .csv file will have two distinct columns and header names. One column will be for the reports before the export name change and the other column will be for the reports created after the export name was changed.

Similar behavior applies if a new data field (with export name) was added to the template.