



## Legacy Uplink

Last Modified on 2022-08-15 18:36

Tricefy Uplink 4.0 has replaced previous versions of Uplink. We will continue to support previous (Legacy) versions of Uplink. Uplink 4.0 is currently only available for Windows.

Legacy Uplink: [https://downloads.tricefy.com/tricefy\\_uplink/trice4.exe](https://downloads.tricefy.com/tricefy_uplink/trice4.exe)

## Downloading & Installing Tricefy Uplink

Unlike Uplink 4.0, there are different Uplinks available depending on your clinic's needs.

### Typical Uplink

Suitable for most small and medium sized clinics - DICOM is anonymized before it's sent to the patient. It is simple and user friendly; no IT resources are needed

 [Typical Uplink Installation Instructions](#)

### Advanced Uplink

Typical Uplink with additional features for integrating with a PACS or EMR:

- Sending studies to Tricefy using ViewPoint
- Routing studies from Tricefy to an external server
- Query and retrieve using a local server

 [Advanced Uplink Instructions](#)

If you are performing clinical studies or have a need for additional anonymization, contact our Customer Success Team (by selecting the [Contact](#) button below) to inquire about our Anonymized Uplink. This Uplink anonymizes data before uploading to Tricefy and discards any non-DICOM content.

## Learn More about Tricefy Uplink:

 [Tricefy Uplink Requirements](#)

 [Making sure your IP address is Static](#)

 [Instructions for a New Computer](#)

 [Instructions for using Uplink to Download from Tricefy](#)

 [Learn more about anonymization](#)



Did you ever wonder if you could install Uplink on your Android phone?  
[There is an app for that!](#)

Instead of uploading exams to Tricefy, you can also download exams using the Tricefy Downlink



 [Instructions for installing Downlink](#)

---