



## Legacy Uplink & Downlink

Last Modified on 2024-02-07 10:55

Tricefy Uplink version 3 has been replaced by [Uplink 4](#). Only install this older version if needed.



### Legacy Uplink (v3) Download

• [https://downloads.tricefy.com/tricefy\\_uplink/trice4.exe](https://downloads.tricefy.com/tricefy_uplink/trice4.exe)

## Installation Instructions

There are different installation options available with legacy uplink depending on your needs.

### [Typical Uplink](#)

Suitable for most small and medium sized clinics - DICOM is anonymized before it's sent to the patient. It is simple and user friendly; no IT resources are needed

 [Typical Uplink Installation Instructions](#)

### [Advanced Uplink](#)

Typical Uplink with additional features for integrating with a PACS or EMR:

- Sending studies to Tricefy using ViewPoint
- Routing studies from Tricefy to an external server
- Query and retrieve using a local server

 [Advanced Uplink Instructions](#)

If you are performing clinical studies or have a need for additional anonymization, contact our Customer Success Team (by selecting the [Contact](#) button below) to inquire about our Anonymized Uplink. This Uplink anonymizes data before uploading to Tricefy.

## Learn More about Tricefy Uplink:



[Tricefy Uplink Requirements](#)



[Instructions for using Uplink to Download from Tricefy](#)



[Learn more about anonymization](#)



Did you ever wonder if you could install Uplink on your Android phone?  
There is an app for that!

Instead of uploading exams to Tricefy, you can also download exams using the Tricefy Downlink



 [Instructions for installing Downlink](#)

---