


## Study List Export

Last Modified on 2026-07-08 00:43

You can generate a report (CSV) of the Tricefy study list and other details. This feature allows you report on select fields from study data in order to find some key metrics from the clinical operations.

Note that this feature only can pull from actively available study data in your Tricefy accounts. If you are not using Tricefy for long term storage of your studies (e.g. 90 days), data will only be available for the studies currently in your account.

To generate a study list export go to the *User Menu > Report Exports > Create Study List Export*



In the prompt, select the account(s), date range (based on study date), and fields you want to include in the study list export. The report will always include Tricefy Account ID and Account Name for in each row as the first 2 columns/fields.

### Create New Study List Report

Account: Trice Demo

Study date range: Start date 07/01/2026, End date 07/31/2026

Select fields to include in study list report

Study List Fields	Other Study Details
<input checked="" type="checkbox"/> Patient ID	<input checked="" type="checkbox"/> Tricefy Study ID
<input type="checkbox"/> Full Name	<input type="checkbox"/> DICOM StudyInstanceUID
<input checked="" type="checkbox"/> Study Date	<input checked="" type="checkbox"/> Count of study instances/files
<input type="checkbox"/> Patient DOB	<input checked="" type="checkbox"/> Count of patient shares
<input type="checkbox"/> Gender	<input type="checkbox"/> Count of Collaborations
<input type="checkbox"/> Procedure	<input type="checkbox"/> Finalized Report
<input type="checkbox"/> Labels	<input type="checkbox"/> Finalized Report Name
<input type="checkbox"/> Study Created	<input type="checkbox"/> Finalized Report Datetime
<input type="checkbox"/> Study Updated	<input type="checkbox"/> Finalized by User
<input type="checkbox"/> Status	
<input type="checkbox"/> Accession Number	
<input type="checkbox"/> Source	
<input type="checkbox"/> Destination	
<input type="checkbox"/> Sonographers	
<input type="checkbox"/> Ref. Phys	
<input type="checkbox"/> Perf. Phys	
<input type="checkbox"/> Volume Data	
<input type="checkbox"/> Interface Status	

Cancel Begin Export

After clicking "Begin Export", the request will be added in the export list. Once the report is ready a downloadable link will be available as a .zip file. Unpacking the zip file will contain the CSV which you can import into any spreadsheet tool for further analysis and filtering.